Scope and Access


The Sequence Viewer (SV) from NCBI is a free, web-based tool for graphical display and analysis of genomic, mRNA, and protein sequences. It is a flexible interface that combines a bird's eye view of a complete chromosome with a detailed interrogation at a single nucleotide resolution. You can access this display from individual records in Nucleotide, Protein, Gene, Genome, SNP, dbGaP, dbVar and Clone databases, and from various tools such as BLAST and Primer BLAST. You can also import and render custom data in different formats into SV. For more details see the SV help document linked on the project's homepage at https://www.ncbi.nlm.nih.gov/projects/sviewer/



Contact: info@ncbi.nlm.nih.gov


Activating an SV Display Through the *Graphics* Link

A common use of SV is to display a large genome. For example, you can retrieve the human chromosome 1 record from the Nucleotide database by its accession number NC_000001, and display it graphically in SV by selecting <u>Graphics</u> under the display format menu under <u>GenBank</u> (A) or simply by clicking the <u>Graphics</u> link (B) at the top of the record.




SV Display: the Overview Panel

The top-most panel of the Graphics page represents chromosome 1 in its entirety (**C**) with its ruler showing the genomic coordinates. You can search for Genes and other annotated features using the "Search" option in the "Tools" dropdown menu (**D**), where the example searches for the IL6R (**E**) gene. Click the "IL6R" gene (**F**) to zoom in to that feature for more detailed examination (p. 2).




SV Display: the Graphical Panel

The graphical panel below the overview shows the details of the IL6R gene, with its genomic range shown in the ruler at the top of the panel (A). Click the green Gene bar to see details in the expanded section, which depicts exons as vertical lines/boxes (B) linked by introns (as thin lines, C). Arrows (D) in the display indicates the gene orientation. Hover over a displayed item, such as a specific transcript variant (E), to see additional details in a popup.




(F) for manipulating the display. Use the "Set New Marker At Position" option (G) to add a marker, such as "TStop_Variant3," at a desired location. The marker context menu enables single-click "Zoom To Sequence At Marker" function (H). Use the zoom slider (I) to zoom in and out incrementally. Pan the display to surrounding regions by holding the left mouse button and dragging the image horizontally, or by using the arrows at the top (J). Click the "Tracks" button (K) and select "Configure Tracks" to access available tracks in the "Configure Page" dialog box (L), through categories listed under the "Tracks" tab, or locate specific ones using the "Search tracks" function (M). Embedded SV displays from records in other NCBI database, such as records from the Gene database, provides NCBI recommended track sets in a cascading menu under Tracks (K). You can also save your track setting to My NCBI for reference and easy access.




Examining a Specific Region at Single Nucleotide Resolution

A marker place at a position of interest (A) enables ready access to sequences through the "Zoom to Sequence At Marker" option in the marker context menu (B). The zoomed in display shows the CDS stop for the transcript variant 3 of the



Additional Displays and Functions

Using the "Tracks" dialog box (A), you can also upload custom data from external sources into a Sequence Viewer display, such as a BLAST result using the assigned Request ID (B) or data file from other sources in various formats (C).



http://bit.ly/ncbi_factsheets