



# Using SRA Toolkit to BLAST SRA Data Locally

Using blastn\_vdb and tblastn\_vdb to search SRA data locally

<https://www.ncbi.nlm.nih.gov/Traces/sra/?view=software>

National Center for Biotechnology Information • National Library of Medicine • National Institutes of Health • Department of Health and Human Services

## Overview

The ongoing revolution in DNA sequencing technology has produced huge quantities of sequencing data (next-gen sequence data) deposited to public databases. At NCBI, next-gen sequence data are available from the Sequence Read Archive (SRA, 1) database. The SRA toolkit, a suite of clients and standalone tools for working with SRA data, provides convenient access to these sets of data. This toolkit includes BLAST programs (blastn\_vdb and tblastn\_vdb) for aligning sequences against next-gen sequence data. These BLAST programs search SRA formatted databases (VDB) directly and can work as clients to access data at NCBI or downloaded and stored on your local disk. This handout demonstrates how to download and set-up the SRA toolkit on a PC running Windows 7, with the goal of using the SRA BLAST programs to search next-gen data, plus the contigs from whole genome shotgun (WGS) and transcriptome shotgun assemblies (TSA), which are also stored in the SRA native format here at NCBI.

## Downloading and Installing the SRA Toolkit

You can access the download page from the NCBI SRA homepage by clicking the “Download SRA Toolkit” link (A). On the download page, right-click the “MS Windows 64 bit architecture” link (B), and select “Save link as ...” option (C) to save the package to a desired location with enough disk space to accommodate SRA datasets needed. In this example, the package is saved to the desktop (D) using the default name.

The screenshot illustrates the steps to download and install the SRA Toolkit. It shows the SRA homepage, the download page, the selection of the Windows 64-bit architecture, the saving of the toolkit as a zip file to the desktop, and the unzipping of the file.

Right click the icon of the downloaded zip file (E), hover over the “WinZip” option, and select “Unzip to here” option (F) to extract and install the package (G) to Desktop. This creates the sratoolkit.###-win64 folder (directory) on the desktop.

The ### is the version number of the package. It is 2.9.0 at the preparation of this handout. Your installation will be of a newer version.

## Configuring the SRA Toolkit

This configuration allows sratoolkit programs to store downloaded SRA data files in a structured directories and manage remote access to NCBI. This is particularly important for authorized access to encrypted data from dbGaP, since encryption key (.ngc) file and encrypted data are stored and accessed from a defined directory as specified. Steps are:

- Launch the Command terminal by searching with "cmd" and clicking cmd.exe (A).
- In the terminal window, change working directory to the bin subdirectory under the SRA Toolkit (B).
- Type "vdb-config -i" (C) to launch the configuration dialog box.
- Use tab key or specific number key (D) to navigate among the fields. The active field is marked by red (E).
- Press the "Enter" key to toggle the setting on or off for the selected field, with "on" setting marked by "X" (F).
- Tab to the "Change" field and hit the "Enter" key to bring up the "select directory" dialog box (G) to see/adjust the current Public directory structure where SRA, WGS, or RefSeq files will be stored, respectively.
- Press "6" to save the changes, and press "7" to exit and return to terminal window.

The displayed above configuration enables remote access to retrieve data automatically from NCBI, and sets local storage of downloaded file the directory specified by the Public field, "C:\Users\tao\ncbi\public" in this case. To avoid accidental deletion of programs or mixing them with input files, do not use the bin directory as working directory.

Do the following to set the path variable to allow for invoking toolkit programs without directory prefix:

- Click the windows icon at lower left
- Search with environment and select "Edit environment variables for your account" (H)
- Select "path" in "User variable ..." section, click "Edit ..." button (I)
- In the prompted dialog box, click "Edit text ..." button (J)
- Terminate the existing value with a semicolon, then append the path to sratoolkit bin directory and click "OK" (K)

## Example BLAST Searches

The following examples use the example installation shown in p. 1 and 2, with sratoolkit-2.9.0-win64 as the working directory. To do so, launch the command prompt (Start >> Run ... >> typing in "cmd" >> clicking "OK") and change the working directory (cd "PATH\sratoolkit-2.9.0-win64"). Since the absolute directory path varies for installations, the custom portion "C:\users\tao\desktop" of the example command lines is replaced with **PATH** to reduce confusion.

### Example One: Surveying the abundance of *Prochlorococcus marinus* at different depths of the ocean

**Background:** In the ocean, the intensity of available light decreases as depth increases. The density of primary producers decreases due to reduced light available for photosynthesis. The cyanobacterium *Prochlorococcus marinus* is the most abundant primary producer in the open ocean [a] with full genome available [d]. An ocean metagenome project is also available, which systematically sequenced biosamples collected from different depths [c]. This example, adopted from an NCBI news entry ([https://www.ncbi.nlm.nih.gov/books/NBK431007/#news\\_11-19-2013-SRA-BLAST](https://www.ncbi.nlm.nih.gov/books/NBK431007/#news_11-19-2013-SRA-BLAST)), uses the PsaAB sequence from *Prochlorococcus marinus* as query and the blastn\_vdb program to quantify the relative abundance of these genes in sequence reads from different depths.

**Search Setup:** Right click [AF180967](#) and save it to the sratoolkit-2.9.0-win64 directory for use as the query. Table 1 lists a set of SRR accessions corresponding to reads from different depths for use as blast databases. Run the command below to do the search (A). From left to right, the command switches instruct the machine to: run blastn\_vdb program in discontinuous megablast mode, use P\_marinus\_Psa.nt as the query, search against the specified datasets (red), use no dust filter, limit the significance of saved hits to Expect value of 0.1 or better, set the database size to a fixed value (for all searches), asks for maximal 5000 hits (if there are that many) in tabular output, and save the results in designated file (red). Replace input for db and out to do the rest.

**Result Checking:** Using expect value below  $1 \times 10^{-3}$  (reported as #e-#, with # being any digit) as a rough estimate, we can approximate significant hits from each result

```
blastn_vdb -task dc-megablast -query P_marinus_Psa.nt -db
"SRR020493 SRR020494" -dust no -evalue 0.1 -dbsize 300000000
-max_target_seqs 5000 -outfmt 6 -out
P_marinus_enzyme_025m.tab
```

```
find /C "e-" P_marinus_enzyme_025m.tab
----- P_MARINUS_ENZYME_025M.TAB: 362
```

file to see the relative abundances of this organism (B).

From the summary in Table 2, we can see that the abundance of these photosynthetic genes peaks at 75m depth and drops significantly at 500m.

### Example Two: The level of expression of ftsA gene in *Lactococcus piscium* at different growth points

**Background:** *L. piscium* is an organism involved in food spoilage and has been the focus of several published studies including the one with data deposited in SRA (Table 3). We will use the abundance of a cell division protein, ftsA, as an indicator of active cell division and growth.

**Search Setup:** Right-click [CEN28187.1](#) to save it to the sratoolkit-2.9.0-win64 directory. Run the search using the command below (C), which instructs the computer to: run tblastn\_vdb program, use L\_piscium\_fstK.aa as the query to search against datasets in quotes (red), use a stringent word size of 6 without SEG filter, save only hits with significance (Expect value) of 0.1 or better in a database size of 4 billion, get maximal 10000 hits (if there are that many) in tabular format, and save the results in the named file (red).

**Result Checking:** We use the same approach described in Example One to count the significant matches from each time point (D). From the summary in Table 4, we can see that significant hits for fstK gene peaks at Hour 5 and drops significantly by Hour 11. The finding is consistent with the need of exponential growth phase, when ftsK needs to be expressed at a higher level in preparation for cell division.

```
tblastn_vdb -query L_piscium_fstK.aa -db "ERR739203 ERR739201 ERR739206" -word_size 6 -seg no
-evalue 0.1 -dbsize 4000000000 -max_target_seqs 20000 -outfmt 6 -out L_piscium_fstK_3hr.tab
```

```
find /C "e-" L_piscium_fstK_3hr.tab
----- L_PISCIMUM_FSTK_3HR.TAB: 10254
```

Table 1. Available data from an ocean depth study ([PRJNA16339](#))

| Experiment Accession & Title | Run Accession (as db) | Spots   |
|------------------------------|-----------------------|---------|
| SRX007372, HOT186_25m_gDNA   | SRR020493 SRR020494   | 623,559 |
| SRX007369, HOT186_75m_gDNA   | SRR020488 SRR020489   | 673,674 |
| SRX007370, HOT186_110m_gDNA  | SRR020490             | 473,116 |
| SRX007371, HOT186_500m_gDNA  | SRR020491 SRR020492   | 995,747 |

Table 2. Result summary for Example One

| Depth | Significant matches | Spots   | Matches/spots (%) |
|-------|---------------------|---------|-------------------|
| 25m   | 362                 | 623,559 | 0.058             |
| 75m   | 460                 | 673,674 | 0.068             |
| 110m  | 158                 | 473,116 | 0.034             |
| 500m  | 42                  | 995,747 | 0.004             |

Table 3. *L. piscium* Time Course Data ([PRJEB8313](#))

| Experiment Accessions | Time Point | Run (Read) Accessions | Spots      |
|-----------------------|------------|-----------------------|------------|
| ERX682888             | 3 hr       | ERR739203             | 38,698,870 |
| ERX682890             |            | ERR739201             |            |
| ERX682882             |            | ERR739206             |            |
| ERX682883             | 5 hr       | ERR739205             | 4,306,066  |
| ERX682887             |            | ERR739207             |            |
| ERX682889             |            | ERR739199             |            |
| ERX682884             | 11 hr      | ERR739202             | 46,691,504 |
| ERR682886             |            | ERR739200             |            |
| ERR682885             |            | ERR739204             |            |

Table 4. Result summary for Example Two

| Time Point | Significant Matches | Spots      | Matches/Spots (%) |
|------------|---------------------|------------|-------------------|
| 3 hr       | 11297               | 38,698,870 | 0.000292          |
| 5 hr       | 11764               | 4,306,066  | 0.000273          |
| 11 hr      | 6340                | 46,691,504 | 0.000136          |



## Example BLAST Searches (cont.)

In addition to SRA, whole genome shotgun contigs (WGS) and transcriptome shotgun sequence assembly (TSA) also adopt the VDB format to store sequences. This enables access to these datasets through tools provided by the *sratoolkit*. When a specific dataset needs to be accessed repeatedly for a series of analyses, downloading and archiving the dataset locally save band-width and make searches go faster.

### Example Three: Finding the MLH1 homolog from *Hydra vulgaris* using downloaded WGS & TSA datasets

**Background:** *H. vulgaris* is an important model organism in biological research. Unfortunately, genomic and transcript sequences available for this organism are often partial or lack annotation.

**Search Setup:** This example finds a full-length transcript from TSA data and the genomic sequence from WGS.

- **Save the query:** Right-click the MLH1 protein sequence

- **NP\_499796** from *C. elegans* to save it to the *sratoolkit* directory

- **Download the database:** Run *prefetch* to download the datasets for projects GAOL (TSA) and ACZU (WGS) (A), which takes advantage of the *fasp* protocol (blue, Aspera installation required) and saves the file to the "PATH\ncbi\public\wgs\" directory. Summaries of these two project are available online at: <http://1.usa.gov/1b6FBn9>.

- **Run *tblastn\_vdb*:** Search the TSA dataset with *tblastn\_vdb* to locate the assembled transcript for this gene from *Hydra magnipapillata* and saves the tabular formatted output to a file (B)

- **Examine the output:** Page through the result with "more" command, the first match (blue) is most likely the transcript for MLH1 homolog from *Hydra magnipapillata*.

- **Retrieve the top match,** [GAOL01023314.1](#) (right click and select "save link as ...") to save it to SRA Toolkit directory, then search it against the WGS dataset to locate the genomic counterpart(s) using *blastn\_vdb* (C).

**Result Checking:** The final *blastn\_vdb* search results (D) indicate that the WGS assembly is incomplete for the MLH1 transcript matched to two contigs, the 5'-portion matched to ACZU01088177 (red) and re-

maining 3'- portion matched to ACZU010338771 (blue), respectively. This piece of information could be used to place these contigs into a longer scaffold. These matches also reveal the general exonic structure of this gene. Note that we can use "-outfmt 7" instead of "-outfmt 6" to see the column header description.

```
PATH\sratoolkit.2.9.0-win64>prefetch GAOL01
Maximum file size download limit is 20,971,520KB

2018-05-29T17:57:13 prefetch.2.9.0: 1) Downloading 'GAOL01'...
2018-05-29T17:57:13 prefetch.2.9.0: Downloading via fasp...
2018-05-29T17:57:17 prefetch.2.9.0: fasp download succeed
2018-05-29T17:57:17 prefetch.2.9.0: 1) 'GAOL01' was downloaded successfully
2018-05-29T17:57:17 prefetch.2.9.0: 'GAOL01' has 0 unresolved dependencies
```

```
PATH\sratoolkit.2.9.0-win64>prefetch ACZU01
Maximum file size download limit is 20,971,520KB

2018-05-29T17:57:17 prefetch.2.9.0: 2) Downloading 'ACZU01'...
2018-05-29T17:57:17 prefetch.2.9.0: Downloading via fasp...
2018-05-29T17:57:24 prefetch.2.9.0: fasp download succeed
2018-05-29T17:57:24 prefetch.2.9.0: 2) 'ACZU01' was downloaded successfully
2018-05-29T17:57:24 prefetch.2.9.0: 'ACZU01' has 0 unresolved dependencies
```

```
PATH\sratoolkit.2.9.0-win64>dir C:\users\tao\ncbi\public\wgs
[... ..]
Directory of PATH\ncbi\public\wgs
```

```
Directory of C:\users\tao\ncbi\public\wgs

05/29/2018 01:57 PM <DIR> .
05/29/2018 01:57 PM <DIR> ..
11/18/2017 12:28 PM 210,134,171 ACZU01
02/17/2017 04:24 AM 18,077,882 GAOL01
                2 File(s) 228,212,053 bytes
                2 Dir(s) 358,024,933,376 bytes free
```

```
PATH\sratoolkit.2.9.0-win64>tblastn_vdb -query C_elegans_MLH1.aa -db
GAOL01 -seg no -evaluate 0.01 -max_target_seqs 1000 -outfmt 6 -out
Hydra_MLH1.tab
```

```
PATH\sratoolkit.2.9.0-win64>more Hydra_MLH1.tab
gi|71991825|ref|NP_499796.2| gi|550390740|gb|GAOL01023314.1| 34.16 767
425 23 4 756 54 2156 7e-114 364
gi|71991825|ref|NP_499796.2| gi|550392993|gb|GAOL01022348.1| 30.40 352
213 12 4 330 216 1250 3e-033 137
gi|71991825|ref|NP_499796.2| gi|550381927|gb|GAOL01025681.1| 26.56 384
265 8 2 375 74 1204 2e-031 131
```

```
PATH\sratoolkit.2.9.0-win64>blastn_vdb -query Hydra_MLH1_transcript.nt
-db ACZU01 -dust no -evaluate 0.01 -outfmt 6 -out Hydra_MLH1_contigs.tab
```

```
PATH\sratoolkit.2.9.0-win64>more Hydra_MLH1_contigs.tab
gi|550390740|gb|GAOL01023314.1| gi|261986949|gb|ACZU01033877.1| 100.00 650
0 0 488 1137 18087 17438 0.0 1201
gi|550390740|gb|GAOL01023314.1| gi|261986949|gb|ACZU01033877.1| 99.45 548
2 1 1360 1906 5410 4863 0.0 994
gi|550390740|gb|GAOL01023314.1| gi|261986949|gb|ACZU01033877.1| 99.27 273
2 0 1900 2172 2509 2237 3e-137 494
gi|550390740|gb|GAOL01023314.1| gi|261986949|gb|ACZU01033877.1| 99.46 184
1 0 1138 1321 6671 6488 2e-089 335
gi|550390740|gb|GAOL01023314.1| gi|261986949|gb|ACZU01033877.1| 100.00 40
0 0 1320 1359 5539 5500 4e-011 75.0
gi|550390740|gb|GAOL01023314.1| gi|261932540|gb|ACZU01088177.1| 100.00 478
0 0 15 492 2281 1804 0.0 883
```

## Technical Assistance

Send BLAST-related comments, questions, and bug reports to: [blast-help@ncbi.nlm.nih.gov](mailto:blast-help@ncbi.nlm.nih.gov)

Send other non-BLAST related questions to: [info@ncbi.nlm.nih.gov](mailto:info@ncbi.nlm.nih.gov)